



GENO-N

GENOLUTION NGS LIBRARY PREP



Automated library preparation by **Geno-N**





Genolution, a leading molecular diagnostics company, is launching advanced NGS Library Preparation Equipment. This technology fully automates the preparation of fragmented DNA for Next Generation Sequencing (NGS).



FEATURES



- 1 Processes up to 96 sample libraries at once
- 2 Reduces labor time and boosts diagnostic accuracy
- 3 Minimizes human error
- 4 Reduces reagent loss
- 5 Enhances inspection convenience with user-friendly operation

SPECIFICATION

	Component	Description
8-ch Pipetting module 	8-ch Pipetting module	CV% Limit 10 ul : <2% (10 ul tip) 50 ul : <0.5% (300 ul tip) 100 ul : <0.5% (300 ul tip)
Heating & Shaking module 	Heating & Shaking module	Temp.(°C) : RT~100 RPM : 0~1600
Flexible Reservoir rack 	Flexible Reservoir rack	Top plate can be manufactured according to user's purpose
Tip types 	Tip types	300 ul, 50 ul, 10ul
Dimension (W * D * H) 	Dimension (W * D * H)	900 * 750 * 850 (mm)
Weight 	Weight	85 kg